
Avalon Elite

Bi-Caval Dual Lumen Catheter







AVALONELITE®

The Avalon Elite advantages:

- Inserted into the internal jugular vein this patented device is able to match the body's natural flow ratios by simultaneously removing blood from both, the SVC and IVC and returning blood to the RA
- Large family of sizes broadens clinical applications for neonatal, pediatric, and adult patients
- Constructed with an exclusive material that combines the durability of polyurethane and the flexibility and biostability of silicone
- Radiopaque to assist in catheter insertion and placement

Featuring the Avalon Elite introducer system

- "Fingerlike" catheter tip reinforcements for optimal catheter introducer interface
- Tapered introducer tip is designed for improved percutaneous insertion
- Ability to customize introducer tip extension based on physician preference



Reinfusion port



Catheter tip with introducer


Avalon Elite Bi-Caval Dual Lumen Catheters

Large family of sizes for broad clinical application.

Pressure drop vs. flow 13, 16, 19 Fr.

250 19 Fr. 200 150 100 Pressure drop [mmHg] -100 -150 2.0 Flow Rate I/min (H₂O at ambient temperature)

Pressure drop vs. flow 20, 23, 27, 31 Fr.



Avalon Elite Bi-Caval Dual Lumen Catheters

| Article No. | Size | Insertable Length | Connector Size | Order Unit |
|-------------|------------------|-------------------|----------------|------------|
| 10013 | 13 Fr. (4.3 mm) | 11 cm (4.3") | 1/4" | 1/Carton |
| 10016 | 16 Fr. (5.3 mm) | 14 cm (5.5") | 1/4" | 1/Carton |
| 10019 | 19 Fr. (6.4 mm) | 21 cm (8.3") | 1/4" | 1/Carton |
| 10020 | 20 Fr. (6.7 mm) | 31 cm (12.2") | 3/8" | 1/Carton |
| 10023 | 23 Fr. (7.7 mm) | 31 cm (12.2") | 3/8" | 1/Carton |
| 10027 | 27 Fr. (9.0 mm) | 31 cm (12.2") | 3/8" | 1/Carton |
| 10031 | 31 Fr. (10.3 mm) | 31 cm (12.2") | 3/8" | 1/Carton |

The Avalon Elite vascular access kits*

for use with catheters and cannulae

The Avalon Elite advantages:

- Extended dilator tip is designed for smooth vessel insertion and guide wire tracking
- Large range of dilator sizes accept standard 0.038" guide wires
- Non-coated stainless steel guide wire displays 10 cm incremental depth markings
- 3 mm J-tip guide wire end can be temporarily straightened to aid placement



100 cm Guide Wire Kit

- 10 Fr., 12 Fr., 16 Fr., 20 Fr., 24 Fr. dilators
- 0.038" x 100 cm guide wire with 10 cm depth mark increments
- 18 ga. puncture needle
- Scalpel blade
- 10 ml syringe

210 cm Guide Wire Kit

- 10 Fr., 12 Fr., 16 Fr., 20 Fr., 24 Fr., 26 Fr., 30 Fr. dilators
- 0.038" x 210 cm guide wire with 10 cm depth mark increments
- 18 ga. puncture needle
- Scalpel blade
- 10 ml syringe

Vascular Access Kits

| Article No. | Guide Wire Length | Order Unit |
|-------------|-------------------|------------|
| 12100 | 100 cm | 5/Carton |
| 12210 | 210 cm | 5/Carton |



This brochure contains information about products which may be pending regulatory approval to be marketed in your country.

Contact your local Maquet representative for more information.

See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

MAQUET GETINGE GROUP

Maquet Cardiopulmonary GmbH

Kehler Str. 31

76437 Rastatt, Germany

Phone: +49 7222 932-0

Fax: +49 7222 932-1888 info.cp@maquet.com

www.maquet.com

GETINGE GROUP

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, Getinge and Maquet. ArjoHuntleigh focuses on patient mobility and wound management solutions. Getinge provides solutions for infection control within healthcare and contamination prevention within life sciences. Maquet specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.